Choix de la langue / Choice of language / Elección de idioma / Taalkeuze / Sprachauswahl

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[→ EU User Guide in English](#)

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
CE marking for European compliance. The application complies with the European Union Directive on medical devices.

## Table of content

<b>1. OdySight presentation .....</b>	<b>3</b>
<b>1.1 Intended use .....</b>	<b>3</b>
<b>1.2 Product presentation .....</b>	<b>3</b>
<b>1.3 Contre-indications.....</b>	<b>4</b>
<b>1.4 Product description .....</b>	<b>4</b>
<b>1.5 Product performance .....</b>	<b>5</b>
<b>1.6 Security .....</b>	<b>7</b>
<b>2. General set up and use .....</b>	<b>8</b>
<b>2.1 Installing OdySight.....</b>	<b>8</b>
<b>2.2 Connect to OdySight .....</b>	<b>8</b>
<b>2.3 Menu and navigation .....</b>	<b>9</b>
<b>2.4 To start a test.....</b>	<b>9</b>
<b>2.5 Test results .....</b>	<b>11</b>
<b>2.6 Close the app .....</b>	<b>11</b>
<b>2.7 Application updates.....</b>	<b>12</b>
<b>3. Technical informations .....</b>	<b>13</b>
<b>3.1 Minimal requirements and compatible devices.....</b>	<b>13</b>
<b>3.2 Battery .....</b>	<b>13</b>
<b>3.3 Screen cleanliness .....</b>	<b>13</b>
<b>4. Warnings and precautions .....</b>	<b>13</b>
<b>5. Contact us .....</b>	<b>14</b>
<b>6. Symbols used on the device labelling and user guide .....</b>	<b>15</b>

## 1. OdySight presentation

---

### 1.1 Intended use

OdySight® is a medical software intended for remote monitoring of the visual parameters: near visual acuity and symptoms such as metamorphopsia and scotoma. The test results of user's self-evaluation are transmitted to the physician.

**OdySight® doesn't interpret the test results and it is not intended to diagnose but rather to help patients determine if they need to consult a physician for diagnosis and treatment.**

Indications: OdySight® is indicated for adult patients aged 18 and over.

### 1.2 Product presentation

OdySight consists of an app for the user and a dashboard for the physician and Healthcare Professionals.

OdySight could be recommended for users that, from the point of view of the ophthalmologist or other healthcare professionals, the user will benefit from remote monitoring of visual parameters: measure of near visual acuity or the evaluation of metamorphopsia and scotomas.

OdySight could be recommended for patients requiring regular, close monitoring of vision. It is to assist the physicians with the patient's visual acuity follow up care between scheduled clinic visits.

The app is prescribed to the patient by a healthcare professional.

OdySight doesn't interpret the test results and is not intended to diagnose. The diagnosis and patient treatment remain the responsibility of the physician.

OdySight does not replace scheduled clinic visits.

**The physician prescribing the device must ensure that the patient is competent to use the device.**

OdySight shall be prescribed to adult users who have:

- a) Binocular visual acuity with habitual correction equal to or better than 3/10 (0,5 logMAR; (20/63))
- b) Access to compatible smartphone or tablet
- c) A valid email address and phone number to create an account
- d) The competence to use the app and complete the test modules remotely, without the guidance of a healthcare professional
- e) The capability to regularly perform a simple self-test at home

### 1.3 Contre-indications

OdySight app is contra-indicated for:

- epileptic users
- binocular visual acuity with habitual correction less than 3/10 (0,5 logMAR; (20/63))

### 1.4 Product description

OdySight is a digital technology system composed of two different components: the OdySight app for users and OdySight dashboard for healthcare professional.

The OdySight app is divided into the following modules:

- the Near Visual Acuity (NVA) module allows the user to test their visual acuity; based on the ETDRS (Early Treatment for Diabetic Retinopathy Study)
- the Amsler Grid (AG) module evaluates metamorphopsia and scotomas; based on the paper Amsler grid
- the authentication module allows users to securely connect to the OdySight application
- the menu module allows for navigation within the OdySight application
- the game and DLC (Downloadable content) modules: related to the puzzles which are designed for the retention of the users

Note: The game module is not part of the OdySight medical device.


The application is available by prescription only and can be used on supported iOS/Android smartphones that meet the minimum technical requirements (see section 3.1).

- [OdySight app](#)

The app consists of two parts: vision tests (visual acuity, Amsler grid), and a game (puzzles). The vision tests are intended to collect test results similar to those usually obtained in a medical center using standard methods.

The visual acuity test assesses near vision, while the Amsler grid test assesses the evolution of scotomas (black spots) and / or metamorphopsia (deformation).

*Near visual acuity and Amsler Grid tests in the OdySight app*

- **OdySight Dashboard**

The dashboard, intended for use by healthcare professionals, is a secure website for viewing the results of vision tests taken by their patients.

## 1.5 Product performance

OdySight allows the remote monitoring of visual parameters.

OdySight is a tool intended to complement your standard medical monitoring.

**The application does not replace prescribed medical visits or examinations by your physician, which patients must continue to undergo.**

Visual acuity is the most commonly used parameter to assess and monitor chronic eye diseases.

The visual tests available in OdySight have been evaluated and compared with standard tests performed in medical centers during a clinical trial. The clinical study demonstrated that the visual acuity tests and the Amsler Grid provided in OdySight showed a strong concordance with tests performed in medical settings (standard methods used – the ETDRS scale for distance and near visual acuity, and the paper-based Amsler Grid).

Each day, OdySight informs you which eye needs to be tested. On average, an OdySight test takes 5 minutes to complete.

If your test results suggest it, OdySight will recommend that you contact your referring healthcare professional. Their contact details are saved and directly accessible from within the app.


These tests allow you to unlock access to puzzles. Game sessions are relatively short so as not to cause you eye strain.

## AMSLER GRID TEST

The Amsler grid was created by Marc Amsler, a Swiss ophthalmologist, in 1945. It is a self-monitoring tool for maculopathies such as AMD and macula edema.

By staring at the center of the grid with one eye open at a set distance (around 40 cm), a person may see metamorphopsia (wavy lines) or scotoma (black patches) appear, which are symptoms of a change to an optical condition. In this case, monitoring by an ophthalmologist is required.

The Amsler Grid does not send notifications advising you to contact your physician, as may be the case with the visual acuity test. When completing the Amsler grid test on OdySight, if you notice the appearance or worsening of any of these anomalies, please contact your physician.

### 1.6 Security

Your health data can be accessed at any time by your ophthalmologist through the online dashboard, which access is both secure and restricted.

These data, including your medical test results, are stored by a certified health data hosting provider that complies with all applicable legal and regulatory requirements. The transmission of your data from the server to the dashboard is secure and encrypted.

## 2. General set up and use

---

### 2.1 OdySight Installation

The application is only available on smartphones or tablets.

To download OdySight, you will need a device with the minimum technical specification's requirements (cf. section 3.1) and an internet connection.

If you are using an **Apple device**:

1. Go to the App Store
2. Click the search bar
3. Write OdySight and click search
4. Click the **download** button beside the OdySight logo.

If you are using an **Android device**:

1. Go to the Play Store.
2. Click the search bar.
3. Write OdySight and click search.
4. Select the OdySight app.
5. Click the *Download* button.

The OdySight application will then be installed on your device.

Downloading OdySight can also be done by clicking on the link in the email or text message received when registering.

### 2.2 Connect to OdySight

Once the app is installed on your smartphone or tablet, you must log in using two separate codes:

- **The one-time code.** This one-time code is sent to you by email by Tilak Healthcare when your ophthalmologist prescribes the application. It is an 8-digit code that is used only once when you log in to OdySight for the first time.
- **The personal password.** You are responsible for choosing this personal password during your first login to the application. This password will be required each time you log in.

When entering your personal password, you can check the box "Remember me". This way, the app will remember your password, so you won't have to enter it every time you log in.

Since the use of the OdySight app is strictly personal, we recommend that you do not check the Remember me box on a device that may be used by anyone other than yourself.

If you have lost your one-time code or if the app does not accept it, please visit the following page: <https://www.odysight.fr/je-nai-pas-recu-ma-clef-chiffree-ou-elle-ne-fonctionne-pas/>

If you have lost your personal password, please visit the following page: <https://www.odysight.fr/jai-oublie-de-mon-code-personnel/>

## 2.3 Menu and navigation

### The Map

The map is the main menu of the app. From this screen, you can start a vision test, launch a puzzle, change settings, share your health insurance card (for patients monitored in France), and access various pieces of information about OdySight. You will land on this screen as soon as you open the app. You can also be redirected clicking on the information menu icon represented by three horizontal lines at the top right of the screen. To return to the map, click on “Return to the map” at the bottom of the menu.

### Menus

OdySight is composed by two menus accessible from the map:

- The vision test menu: This menu will allow you to start a vision test (Visual acuity or Amsler grid).
- The information menu represented by three horizontal lines at the top right of the screen: This menu will grant you access to different menus: the “Help menu”, “Parameters menu” and “Legal menu”.

- The Help menu: this menu grants you access to frequently asked questions (Questions/ answers), to the user guide and allows you to contact the support service.
- The Parameters menu: This menu grants you access to the different available languages.
- The legal menu: This menu grants you access to the terms of use, Privacy Policy and the « Delete my account » button.


## 2.4 To start a test

You can start a test in two ways:

- When you open the application for the first time each day, a window will appear letting you know you need to take your daily test. You can click accept, then follow the instructions shown on the screen.
- When you click on the Vision test menu, a list of tests will appear. Click on a vision test and follow the instructions shown on the screen.

To take a test in the best possible conditions:

1. Make sure your smartphone or tablet is connected to the internet.
2. Ensure that there is plenty of ambient light, but that none of it is shining directly on your face.
3. Place your device on a stable surface.
4. Make sure your device is in portrait orientation, that is to say vertical, as shown in the picture below.
5. Click on the Vision test menu and choose your vision test.
6. If you wear glasses for nearsightedness, put them on.
7. Make sure you are 40 cm from the screen (plus or minus 5 cm)
8. Follow the instructions that are shown on the screen.

Should you have any issues or see an error message during your test, check the FAQ relating to the correct test conditions: <https://www.odysight.fr/en/je-narrive-pas-a-faire-mon-test/>

**The visual acuity test** is presented in the form of a short and interactive test in which the user must indicate the direction of the letter “E” by swiping their finger in the corresponding direction (top, bottom, right or left).

**The Amsler grid test** is a digitalized version of the paper Amsler grid. By using their finger as a pencil, the user can draw any irregularities which they see appear.

A notification will automatically be sent to your device (smartphone or tablet) to remind you that a visual acuity test is available.

To turn off the notifications, you will find the instructions on this page: <https://www.odysight.fr/en/gerer-les-notifications/>

It is advisable to limit your use of the application to a maximum of 90 minutes at a time.

## 2.5 Test results

The results are shown at the end of each visual test:

- The results for the visual acuity test are shown in VA letters. This score varies from 10 to 85. A variation of a few letters from one test to another is considered as not significant. The results shown are for information purposes only.
- The results of the Amsler grid are presented in the form of an image. If you observe anomalies (wavy lines or dark patch) in the Amsler grid, it is advised to contact your physician.

If the internet connection is interrupted while the test is being carried out, the data will be saved and sent as soon as the connection is restored.

The color of paint splatters appearing in visual acuity test is of purely artistic interest and does not indicate a right or wrong answer.

Visual acuity test results may differ from those obtained during a medical consultation. OdySight measures near visual acuity, not distance visual acuity as may be done during a consultation. Additionally, OdySight may use a different unit system than your ophthalmologist.

If a variation in visual acuity is detected, you will be asked to confirm it with a second test, and a third test if the second one confirms the variation. An SMS will be sent to you after each test result indicating a variation in visual acuity, prompting you to complete a confirmation test. After the third test, you will receive both an SMS and a specific email, and your doctor will be notified if the variation is confirmed. You will then be advised to contact your doctor to schedule an appointment.

There are no notifications of anomalies using the Amsler grid sent by the app. If you see any dots changes, you must contact your physician for visit request.

Please contact your physician for any further questions about Amsler Grid and your test results.

## 2.6 Close the app

The way to close the application depends on your device. Here are several methods for closing the OdySight application.

### **If you have an Apple device**

To close OdySight on an Apple device (iPhone or iPad), follow these instructions:

1. Click twice on the home button (the round button below your screen) quickly. The list of opened apps will appear, with the OdySight window first. If your device model doesn't have the home button, put your finger at the bottom of your screen and slide it from bottom to top to display the list of the opened apps.
2. Slide the OdySight window from top to bottom until it disappears.
3. When the window disappears, the app is closed.

### **If you have an Android™ device**

When you open OdySight, the app is put in full screen. Follow these instructions to leave the full screen mode:

- Put your finger at the bottom of the screen and swipe or slide it up. The triangle, round, and square buttons of your Android device will appear.
- Click the square. The window of the app will reduce in size, and you will see your wallpaper as well as the other opened apps (if there are ones).

The method to close an app is not the same on all Android devices, but here is a list of possibilities to close OdySight:

- Click the cross at the top right side of the OdySight window.
- Click the cross at the top left side of the OdySight window.
- Slide the OdySight window from left to right until it disappears.
- Slide the OdySight window from top to bottom until it disappears.

When the window disappears, the app is closed.

## **2.7 Application updates**

Major application updates happen automatically when your device is connected to the internet. When one of these updates is available, you will see a message appear in OdySight letting you know about the update.

If the update did not occur automatically, a screen prompting you to perform the update manually will appear when the app is launched. To perform the update manually, follow these steps:

1. Click on the "Download" button.
2. You will be redirected to your device's store (App Store for Apple devices; Play Store for Android devices).
3. Tap the "Update" button that appears next to the app's icon. This will initiate the update download.

4. To reopen the app after the update, tap the “Open” button from the store or from your device’s app screen.

### 3. Technical informations

---

#### 3.1 Minimal requirements and compatible devices

Any device meeting the criteria below can download the OdySight app:

- Minimum operating system:
  - Android™ 8
  - iOS® 13
- Screen resolution equal to or higher than: :
  - 280 pixels per inch (PPI) for Android™
  - 260 pixels per inch (PPI) for iOS®
- Screen diagonal greater than 4.7”
- Presence of a front-facing camera (located above the device screen)
- Presence of a brightness sensor (for Android™ devices))

This information is also available on the Help Center page: <https://www.odysight.fr/la-liste-dappareils-supportes-suite-retour-test-qa/>

**If one of your eyes has a visual acuity greater than 6/10 (+0.2 logMAR or 20/32), the resolution of the screen must be at least 320 ppi.**

If you have any questions, don't hesitate to send a message to the following email address: supportFR@odysight.app or call us at +331 89 71 13 36.

#### 3.2 Battery

To take the different tests, the screen of your device must admit sufficient light. This is why the phone should be **properly charged** to avoid it going into “low battery” mode.

If the phone’s battery level is below 15%, the app will no longer allow the user to start any tests. The phone must then be recharged in order to continue.

#### 3.3 Screen cleanliness

For the best test performance, it is important that the device’s screen is intact and free of any deep scratches or dents and is clean.

### 4. Warnings and precautions

---

OdySight does not replace the patients’ healthcare facilities visit and tests prescribed by the healthcare professional as part of monitoring of the vision. Continue follow the visits calendar.

OdySight does not provide an interpretation of test results, diagnose or treatment recommendations. It is not intended to provide a diagnosis. Diagnosis remains the responsibility of the physician.

All medical decisions must be made by a healthcare professional.

If you notice a change in your vision and think there is a need, you need to contact the healthcare professional who prescribed OdySight to you. His contact is available in the application.

The app should be used only by the patient whom the application has been prescribed to ensure the reliability of the data provided to the healthcare professional.

Carefully follow the application instructions to test the correct eye (left or right). Testing both eyes at the same time or testing the wrong eye could lead to inaccurate measurement results being sent to your physician.

Carefully follow the instructions relative to distance to screen and ambient light conditions before starting a medical test.

The test results may differ from those obtained during your medical appointment. They should be considered as relative values.

Test results may differ slightly from one device to another. Please continue to follow the instructions carefully to perform your medical tests in good conditions with your new device.

The OdySight Visual acuity decrease notifications are not intended to have a direct impact or change the clinical care flow for patients suffering from a chronic ocular disease, but it is only intended to help healthcare professionals to improve patient follow-up.

During the tests, the patient must not force their sight. During visual acuity tests, a "I don't know" button will appear after a few seconds without touching the screen. By click it, the size of the "E" will increase. This will not distort test results, feel free to use it.

Do NOT use this product to self-diagnose vision problems.

Do NOT leave the smartphone or tablet within reach of children while performing the tests.

Tilak Healthcare is not liable for data collected by mistake on the device due to misuse or malfunction resulting from abuse, accidents, modifications or negligence.

## 5. Contact us

---

If you experience any issues with OdySight, please refer to the Frequently Asked Questions, which contain answers to common inquiries: <https://www.odysight.fr/centre-daide/>.

Should you still need help with an issue, contact us at the following email address: [assistance@odysight.fr](mailto:assistance@odysight.fr) or call us at +331 89 71 13 36.

**Any serious incident occurring in connection with the device must be brought to the attention of the manufacturer and to the competent authority of the Member State in which the user is established.**

## 6. Symbols used on the device labelling and user guide


---

	CE marking for European compliance. The application complies with the European Union Directive on medical devices.
	Fabricant
	Consult the instructions for use
	Caution/warning
	For your information

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